

Antibodies Product Datasheet

Product Name:	Cripto-3 Antibody
Catalog Number:	C30623
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human TDGF1P3. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	TDGF1P3
Protein Name:	Putative teratocarcinoma-derived growth factor 3
Other Name:	TDGF1P3; CRIPTO3; TDGF2; TDGF3; Putative teratocarcinoma-derived growth factor 3; Cripto-3 growth factor; Epidermal growth factor-like crypto protein CR3; Teratocarcinoma-derived growth factor 1 pseudogene 3
Human Swiss Prot No.:	P51864
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cripto-3 Polyclonal Antibody detects endogenous levels of Cripto-3 protein.
Storage and Stability:	-20°C/1 year