

Antibodies Product Datasheet

Product Name:	Ebi3 Antibody
Catalog Number:	C30627
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human EBI3. AA range:101-150
Concentration:	1 mg/ml
Gene Name:	EBI3
Protein Name:	Interleukin-27 subunit beta
Other Name:	EBI3; IL27B; Interleukin-27 subunit beta; IL-27 subunit beta; IL-27B; Epstein-Barr virus-induced gene 3 protein; EBV-induced gene 3 protein
Human Swiss Prot No.:	Q14213
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Ebi3 Polyclonal Antibody detects endogenous levels of Ebi3 protein.
Storage and Stability:	-20°C/1 year