

Antibodies Product Datasheet

Product Name:	Lefty Antibody
Catalog Number:	C30633
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the C-terminal region of human LEFTY1/LEFTY2. AA range:301-350
Concentration:	1 mg/ml
Gene Name:	LEFTY1/LEFTY2
Protein Name:	Left-right determination factor 1/Left-right determination factor 2
Other Name:	LEFTY1; LEFTB; LEFTYB; Left-right determination factor 1; Left-right determination factor B; Protein lefty-1; Protein lefty-B; LEFTY2; EBAF; LEFTA; LEFTYA; TGFB4; PSEC0024; Left-right determination factor 2; Endometrial bleeding-associated factor; Left-ri
Human Swiss Prot No.:	O75610/O00292
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Lefty Polyclonal Antibody detects endogenous levels of Lefty protein.
Storage and Stability:	-20°C/1 year