

Antibodies Product Datasheet

Product Name:	IL-22 Antibody
Catalog Number:	C30635
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the C-terminal region of human IL22. AA range:121-170
Concentration:	1 mg/ml
Gene Name:	IL22
Protein Name:	Interleukin-22
Other Name:	IL22; ILTIF; ZCYTO18; Interleukin-22; IL-22; Cytokine Zcyto18; IL-10-related T-cell-derived-inducible factor; IL-TIF
Human Swiss Prot No.:	Q9GZX6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IL-22 Polyclonal Antibody detects endogenous levels of IL-22 protein.
Storage and Stability:	-20°C/1 year