

Antibodies Product Datasheet

Product Name:	HRP-3 Antibody
Catalog Number:	C30644
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1/100-1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human HDGFRP3. AA range:131-180
Concentration:	1 mg/ml
Gene Name:	HDGFRP3
Protein Name:	Hepatoma-derived growth factor-related protein 3
Other Name:	HDGFRP3; HDGF2; CGI-142; Hepatoma-derived growth factor-related protein 3; HRP-3; Hepatoma-derived growth factor 2; HDGF-2
Human Swiss Prot No.:	Q9Y3E1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	HRP-3 Polyclonal Antibody detects endogenous levels of HRP-3 protein.
Storage and Stability:	-20°C/1 year