

Antibodies Product Datasheet

Product Name:	IFN- β Antibody
Catalog Number:	C30654
Reactivity:	Human, Rat, Mouse
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	IHC-p 1:50-200, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthetic peptide from human protein at AA range: 91-140
Concentration:	1 mg/ml
Gene Name:	IFNB1 IFB IFNB
Protein Name:	Interferon beta (IFN-beta) (Fibroblast interferon)
Other Name:	Interferon beta (IFN-beta) (Fibroblast interferon)
Human Swiss Prot No.:	P01574
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	The antibody detects endogenous IFN- β
Storage and Stability:	-20°C/1 year