

Antibodies Product Datasheet

Product Name:	IL-26 Antibody
Catalog Number:	C30658
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from Interleukin-26 at AA range: 21-70
Concentration:	1 mg/ml
Gene Name:	IL26
Protein Name:	Interleukin-26
Other Name:	IL26; AK155; Interleukin-26; IL-26; Protein AK155
Human Swiss Prot No.:	Q9NPH9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IL-26 Polyclonal Antibody detects endogenous levels of IL-26 protein.
Storage and Stability:	-20°C/1 year