

Antibodies Product Datasheet

Product Name:	BMP-3A Antibody
Catalog Number:	C30667
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1/100-1/300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the C-terminal region of human BMP3. AA range:421-470
Concentration:	1 mg/ml
Gene Name:	BMP3
Protein Name:	Bone morphogenetic protein 3
Other Name:	BMP3; BMP3A; Bone morphogenetic protein 3; BMP-3; Bone morphogenetic protein 3A; BMP-3A; Osteogenin
Human Swiss Prot No.:	P12645
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	BMP-3A Polyclonal Antibody detects endogenous levels of BMP-3A protein.
Storage and Stability:	-20°C/1 year