

Antibodies Product Datasheet

Product Name:	IFN- ϵ Antibody
Catalog Number:	C30668
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the C-terminal region of human IFNE. AA range:159-208
Concentration:	1 mg/ml
Gene Name:	IFNE
Protein Name:	Interferon epsilon
Other Name:	IFNE; IFNE1; Interferon epsilon; IFN-epsilon; Interferon epsilon-1
Human Swiss Prot No.:	Q86WN2
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IFN- ϵ Polyclonal Antibody detects endogenous levels of IFN- ϵ protein.
Storage and Stability:	-20°C/1 year