

Antibodies Product Datasheet

Product Name:	CTACK Antibody
Catalog Number:	C30669
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-1:300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from C-C motif chemokine 27 at AA range: 51-100
Concentration:	1 mg/ml
Gene Name:	CCL27
Protein Name:	C-C motif chemokine 27
Other Name:	CCL27; ILC; SCYA27; C-C motif chemokine 27; CC chemokine ILC; Cutaneous T-cell-attracting chemokine; CTACK; ESkin; IL-11 R-alpha-locus chemokine; Skinkine; Small-inducible cytokine A27
Human Swiss Prot No.:	Q9Y4X3
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CTACK Polyclonal Antibody detects endogenous levels of CTACK protein.
Storage and Stability:	-20°C/1 year