

Antibodies Product Datasheet

Product Name:	ALR Antibody
Catalog Number:	C30675
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-1:300. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human GFER. AA range:51-100
Concentration:	1 mg/ml
Gene Name:	GFER
Protein Name:	FAD-linked sulfhydryl oxidase ALR
Other Name:	GFER; ALR; HERV1; HPO; FAD-linked sulfhydryl oxidase ALR; Augmenter of liver regeneration; hERV1; Hepatopoietin
Human Swiss Prot No.:	P55789
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ALR Polyclonal Antibody detects endogenous levels of ALR protein.
Storage and Stability:	-20°C/1 year