

Antibodies Product Datasheet

Product Name:	IL-1F10 Antibody
Catalog Number:	C30689
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. IHC-p: 1:100-1:300. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from Interleukin-1 family member 10 at AA range: 101-150
Concentration:	1 mg/ml
Gene Name:	IL1F10
Protein Name:	Interleukin-1 family member 10
Other Name:	IL1F10; FIL1T; IL1HY2; FKSG75; Interleukin-1 family member 10; IL-1F10; FIL1 theta; Interleukin-1 HY2; IL-1HY2; Interleukin-1 theta; IL-1 theta
Human Swiss Prot No.:	Q8WWZ1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IL-1F10 Polyclonal Antibody detects endogenous levels of IL-1F10 protein.
Storage and Stability:	-20°C/1 year