

## Antibodies Product Datasheet

|                               |                                                                                                                                      |
|-------------------------------|--------------------------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | 6CKine Antibody                                                                                                                      |
| <b>Catalog Number:</b>        | C30695                                                                                                                               |
| <b>Reactivity:</b>            | Human, Rat, Mouse                                                                                                                    |
| <b>Applications:</b>          | IHC-p, IF(paraffin section), ELISA                                                                                                   |
| <b>Dilutions:</b>             | IHC-p 1:50-200, ELISA 1:10000-20000                                                                                                  |
| <b>Modification:</b>          | Unmodified/Total                                                                                                                     |
| <b>Source:</b>                | Polyclonal Rabbit                                                                                                                    |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.                |
| <b>Immunogen:</b>             | The antiserum was produced against synthesized peptide derived from the Internal region of human CCL21. AA range:71-120              |
| <b>Concentration:</b>         | 1 mg/ml                                                                                                                              |
| <b>Gene Name:</b>             | CCL21 SCYA21 UNQ784/PRO1600                                                                                                          |
| <b>Protein Name:</b>          | C-C motif chemokine 21 (6CKine) (Beta-chemokine exodus-2) (Secondary lymphoid-tissue chemokine) (SLC) (Small-inducible cytokine A21) |
| <b>Other Name:</b>            | C-C motif chemokine 21 (6CKine) (Beta-chemokine exodus-2) (Secondary lymphoid-tissue chemokine) (SLC) (Small-inducible cytokine A21) |
| <b>Human Swiss Prot No.:</b>  | O00585                                                                                                                               |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                                              |
| <b>Specificity:</b>           | The antibody detects endogenous 6CKine                                                                                               |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                                                         |