

## Antibodies Product Datasheet

|                               |                                                                                                                       |
|-------------------------------|-----------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | HDGF Antibody                                                                                                         |
| <b>Catalog Number:</b>        | C30699                                                                                                                |
| <b>Reactivity:</b>            | Human, Rat, Mouse                                                                                                     |
| <b>Applications:</b>          | IHC-p, IF(paraffin section), ELISA                                                                                    |
| <b>Dilutions:</b>             | IHC-p 1:50-200, ELISA 1:10000-20000                                                                                   |
| <b>Modification:</b>          | Unmodified/Total                                                                                                      |
| <b>Source:</b>                | Polyclonal Rabbit                                                                                                     |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| <b>Immunogen:</b>             | Synthetic peptide from human protein at AA range: 141-190                                                             |
| <b>Concentration:</b>         | 1 mg/ml                                                                                                               |
| <b>Gene Name:</b>             | HDGF HMG1L2                                                                                                           |
| <b>Protein Name:</b>          | Hepatoma-derived growth factor (HDGF) (High mobility group protein 1-like 2) (HMG-1L2)                                |
| <b>Other Name:</b>            | Hepatoma-derived growth factor (HDGF) (High mobility group protein 1-like 2) (HMG-1L2)                                |
| <b>Human Swiss Prot No.:</b>  | P51858                                                                                                                |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                               |
| <b>Specificity:</b>           | The antibody detects endogenous HDGF                                                                                  |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                                          |