

Antibodies Product Datasheet

Product Name:	BMP-6 Antibody
Catalog Number:	C30704
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human BMP6. AA range:341-390
Concentration:	1 mg/ml
Gene Name:	BMP6
Protein Name:	Bone morphogenetic protein 6
Other Name:	BMP6; VGR; Bone morphogenetic protein 6; BMP-6; VG-1-related protein; VG-1-R; VGR-1
Human Swiss Prot No.:	P22004
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	BMP-6 Polyclonal Antibody detects endogenous levels of BMP-6 protein.
Storage and Stability:	-20°C/1 year