

Antibodies Product Datasheet

Product Name:	IL-17B Antibody
Catalog Number:	C30705
Reactivity:	Human, Rat
Applications:	IHC-p, IF(paraffin section), ELISA
Dilutions:	IHC-p 1:50-200, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthetic peptide from human protein at AA range: 1-50
Concentration:	1 mg/ml
Gene Name:	IL17B IL20 NIRF ZCYTO7 UNQ516/PRO1031
Protein Name:	Interleukin-17B (IL-17B) (Cytokine Zcyto7) (Interleukin-20) (IL-20) (Neuronal interleukin-17-related factor)
Other Name:	Interleukin-17B (IL-17B) (Cytokine Zcyto7) (Interleukin-20) (IL-20) (Neuronal interleukin-17-related factor)
Human Swiss Prot No.:	Q9UHF5
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	The antibody detects endogenous IL-17B
Storage and Stability:	-20°C/1 year