

Antibodies Product Datasheet

Product Name:	TLE1/2/3/4 Antibody
Catalog Number:	C30796
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the C-terminal region of human TLE1/TLE2/TLE3/TLE4. AA range:721-770
Concentration:	1 mg/ml
Gene Name:	
Protein Name:	TLE1/2/3/4
Other Name:	similar to transducin-like enhancer of split 1/2/3/4
Human Swiss Prot No.:	Q04724/Q04725/Q04726/Q04727
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	TLE1/2/3/4 Polyclonal Antibody detects endogenous levels of TLE1/2/3/4
Storage and Stability:	-20°C/1 year