

Antibodies Product Datasheet

Product Name:	Calretinin Antibody
Catalog Number:	C30799
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from calbindin 2 at AA range: 191-240
Concentration:	1 mg/ml
Gene Name:	CALB2
Protein Name:	calbindin 2
Other Name:	Calretinin (CR) (29 kDa calbindin)
Human Swiss Prot No.:	P22676
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Calretinin Polyclonal Antibody detects endogenous levels of calbindin 2
Storage and Stability:	-20°C/1 year