

Antibodies Product Datasheet

Product Name:	Ini1 Antibody
Catalog Number:	C30801
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthetic peptide from human protein at AA range: 331-380
Concentration:	1 mg/ml
Gene Name:	SMARCB1 BAF47 INI1 SNF5L1
Protein Name:	SWI/SNF related, matrix associated, actin dependent regulator of chromatin, subfamily b, member 1
Other Name:	SMARCB1 BAF47 INI1 SNF5L1
Human Swiss Prot No.:	Q12824
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	The antibody detects endogenous Ini1 protein
Storage and Stability:	-20°C/1 year