

Antibodies Product Datasheet

Product Name:	Perforin 1 Antibody
Catalog Number:	C30848
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC, ELISA
Dilutions:	WB 1:500-2000, IHC 1:50-200, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the C-terminal region of human PRF1. AA range:451-500
Concentration:	1 mg/ml
Gene Name:	PRF1
Protein Name:	Perforin 1
Other Name:	Perforin-1 (P1) (Cytolysin) (Lymphocyte pore-forming protein) (PFP)
Human Swiss Prot No.:	P14222
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Perforin 1 Polyclonal Antibody detects endogenous levels of Perforin 1
Storage and Stability:	-20°C/1 year