

Antibodies Product Datasheet

Product Name:	IgM Chain C Antibody
Catalog Number:	C30863
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	IHC-p: 100-300.WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from IgM Chain C at AA range: 391-440
Concentration:	1 mg/ml
Gene Name:	IGHM
Protein Name:	IgM Chain C
Other Name:	Ig mu chain C region
Human Swiss Prot No.:	P01871
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IgM Chain C Polyclonal Antibody detects endogenous levels of IgM Chain C
Storage and Stability:	-20°C/1 year