

Antibodies Product Datasheet

Product Name:	CYFIP2 Antibody
Catalog Number:	C30899
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from CYFIP2 at AA range: 1171-1220
Concentration:	1 mg/ml
Gene Name:	CYFIP2
Protein Name:	CYFIP2
Other Name:	Cytoplasmic FMR1-interacting protein 2 (p53-inducible protein 121)
Human Swiss Prot No.:	Q96F07
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CYFIP2 Polyclonal Antibody detects endogenous levels of CYFIP2
Storage and Stability:	-20°C/1 year