

Antibodies Product Datasheet

Product Name:	Integrin β 7 Antibody
Catalog Number:	C30902
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human Integrin β 7. at AA range: 671-720
Concentration:	1 mg/ml
Gene Name:	ITGB7
Protein Name:	Integrin β 7
Other Name:	Integrin beta-7 (Gut homing receptor beta subunit)
Human Swiss Prot No.:	P26010
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Integrin β 7 Polyclonal Antibody detects endogenous levels of Integrin β 7
Storage and Stability:	-20°C/1 year