

Antibodies Product Datasheet

Product Name:	GCK/GLK Antibody
Catalog Number:	C30946
Reactivity:	Human, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human GCK/GLK. at AA range: 221-270
Concentration:	1 mg/ml
Gene Name:	MAP4K2/MAP4K3
Protein Name:	GCK/GLK
Other Name:	Mitogen-activated protein kinase kinase kinase kinase 2 (EC 2.7.11.1) (B lymphocyte serine/threonine-protein kinase) (Germinal center kinase) (GC kinase) (MAPK/ERK kinase kinase kinase 2) (MEK kinase kinase 2) (MEKKK 2) (Rab8-interacting protein) Mitogen-Q12851/Q8IVH8
Human Swiss Prot No.:	
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	GCK/GLK Polyclonal Antibody detects endogenous levels of GCK/GLK
Storage and Stability:	-20°C/1 year