

Antibodies Product Datasheet

Product Name:	Cytochrome b Antibody
Catalog Number:	C30954
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human Cytochrome b. at AA range: 331-380
Concentration:	1 mg/ml
Gene Name:	MT-CYB
Protein Name:	Cytochrome b
Other Name:	Cytochrome b (Complex III subunit 3) (Complex III subunit III) (Cytochrome b-c1 complex subunit 3) (Ubiquinol-cytochrome-c reductase complex cytochrome b subunit)
Human Swiss Prot No.:	P00156
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Cytochrome b Polyclonal Antibody detects endogenous levels of Cytochrome b
Storage and Stability:	-20°C/1 year