

Antibodies Product Datasheet

Product Name:	C8 α Antibody
Catalog Number:	C30956
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human C8 α . at AA range: 201-250
Concentration:	1 mg/ml
Gene Name:	C8A
Protein Name:	C8 α
Other Name:	Complement component C8 alpha chain (Complement component 8 subunit alpha)
Human Swiss Prot No.:	P07357
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	C8 α Polyclonal Antibody detects endogenous levels of C8 α
Storage and Stability:	-20°C/1 year