

Antibodies Product Datasheet

Product Name:	AVPR1A Antibody
Catalog Number:	C30964
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human AVPR1A. at AA range: 11-60
Concentration:	1 mg/ml
Gene Name:	AVPR1A
Protein Name:	AVPR1A
Other Name:	Vasopressin V1a receptor (V1aR) (AVPR V1a) (Antidiuretic hormone receptor 1a) (Vascular/hepatic-type arginine vasopressin receptor)
Human Swiss Prot No.:	P37288
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	AVPR1A Polyclonal Antibody detects endogenous levels of AVPR1A
Storage and Stability:	-20°C/1 year