

Antibodies Product Datasheet

Product Name:	AFG3L2 Antibody
Catalog Number:	C30965
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human AFG3L2. at AA range: 744-793
Concentration:	1 mg/ml
Gene Name:	AFG3L2
Protein Name:	AFG3L2
Other Name:	AFG3-like protein 2 (EC 3.4.24.-) (Paraplegin-like protein)
Human Swiss Prot No.:	Q9Y4W6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	AFG3L2 Polyclonal Antibody detects endogenous levels of AFG3L2
Storage and Stability:	-20°C/1 year