

Antibodies Product Datasheet

Product Name:	FIG4 Antibody
Catalog Number:	C30968
Reactivity:	Human, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human FIG4. at AA range: 341-390
Concentration:	1 mg/ml
Gene Name:	FIG4 KIAA0274 SAC3
Protein Name:	FIG4
Other Name:	Polyphosphoinositide phosphatase (EC 3.1.3.-) (Phosphatidylinositol 3,5-bisphosphate 5-phosphatase) (SAC domain-containing protein 3)
Human Swiss Prot No.:	Q92562
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of FIG4
Storage and Stability:	-20°C/1 year