

Antibodies Product Datasheet

Product Name:	KCTD7 Antibody
Catalog Number:	C30972
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human KCTD7. at AA range: 181-230
Concentration:	1 mg/ml
Gene Name:	KCTD7
Protein Name:	KCTD7
Other Name:	BTB/POZ domain-containing protein KCTD7
Human Swiss Prot No.:	Q96MP8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	This antibody detects endogenous levels of KCTD7
Storage and Stability:	-20°C/1 year