

Antibodies Product Datasheet

Product Name:	IFITM2 Antibody
Catalog Number:	C30976
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human IFITM2. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	IFITM2
Protein Name:	Interferon-induced transmembrane protein 2
Other Name:	IFITM2; Interferon-induced transmembrane protein 2; Dispanin subfamily A member 2c; DSPA2c; Interferon-inducible protein 1-8D
Human Swiss Prot No.:	Q01629
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	IFITM2 Polyclonal Antibody detects endogenous levels of IFITM2 protein.
Storage and Stability:	-20°C/1 year