

Antibodies Product Datasheet

Product Name:	PGLYRP1 Antibody
Catalog Number:	C30977
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/500 - 1/2000. ELISA: 1/10000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human PGLYRP1. AA range:131-180
Concentration:	1 mg/ml
Gene Name:	PGLYRP1
Protein Name:	Peptidoglycan recognition protein 1
Other Name:	PGLYRP1; PGLYRP; PGRP; TNFSF3L; SBBI68; Peptidoglycan recognition protein 1; Peptidoglycan recognition protein short; PGRP-S
Human Swiss Prot No.:	O75594
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	PGLYRP1 Polyclonal Antibody detects endogenous levels of PGLYRP1 protein.
Storage and Stability:	-20°C/1 year