

## Antibodies Product Datasheet

|                               |                                                                                                                       |
|-------------------------------|-----------------------------------------------------------------------------------------------------------------------|
| <b>Product Name:</b>          | AVP Receptor V2 Antibody                                                                                              |
| <b>Catalog Number:</b>        | C30983                                                                                                                |
| <b>Reactivity:</b>            | Human, Mouse                                                                                                          |
| <b>Applications:</b>          | WB, ELISA                                                                                                             |
| <b>Dilutions:</b>             | WB 1:500-2000, ELISA 1:10000-20000                                                                                    |
| <b>Modification:</b>          | Unmodified/Total                                                                                                      |
| <b>Source:</b>                | Polyclonal Rabbit                                                                                                     |
| <b>Purification:</b>          | The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen. |
| <b>Immunogen:</b>             | Synthesized peptide derived from human AVP Receptor V2. at AA range: 1-50                                             |
| <b>Concentration:</b>         | 1 mg/ml                                                                                                               |
| <b>Gene Name:</b>             | AVPR2                                                                                                                 |
| <b>Protein Name:</b>          | AVP Receptor V2                                                                                                       |
| <b>Other Name:</b>            | Vasopressin V2 receptor (V2R) (AVPR V2) (Antidiuretic hormone receptor) (Renal-type arginine vasopressin receptor)    |
| <b>Human Swiss Prot No.:</b>  | P30518                                                                                                                |
| <b>Formulation:</b>           | Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.                                               |
| <b>Specificity:</b>           | AVP Receptor V2 Polyclonal Antibody detects endogenous levels of AVP Receptor V2                                      |
| <b>Storage and Stability:</b> | -20°C/1 year                                                                                                          |