

Antibodies Product Datasheet

Product Name:	FA2H Antibody
Catalog Number:	C30996
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human FA2H. at AA range: 101-150
Concentration:	1 mg/ml
Gene Name:	FA2H
Protein Name:	FA2H
Other Name:	Fatty acid 2-hydroxylase (EC 1.-.-.) (Fatty acid alpha-hydroxylase)
Human Swiss Prot No.:	Q7L5A8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	FA2H Polyclonal Antibody detects endogenous levels of FA2H
Storage and Stability:	-20°C/1 year