

Antibodies Product Datasheet

Product Name:	C8 β Antibody
Catalog Number:	C31006
Reactivity:	Human, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human C8B. AA range:371-420
Concentration:	1 mg/ml
Gene Name:	C8B
Protein Name:	C8 β
Other Name:	Complement component C8 beta chain (Complement component 8 subunit beta)
Human Swiss Prot No.:	P07358
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	C8 β Polyclonal Antibody detects endogenous levels of C8 β
Storage and Stability:	-20°C/1 year