

Antibodies Product Datasheet

Product Name:	FBXO7 Antibody
Catalog Number:	C31009
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthetic peptide from human protein at AA range: 371-420
Concentration:	1 mg/ml
Gene Name:	FBXO7 FBX7
Protein Name:	F-box protein 7
Other Name:	FBXO7 FBX7
Human Swiss Prot No.:	Q9Y3I1
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	The antibody detects endogenous FBXO7 protein
Storage and Stability:	-20°C/1 year