

Antibodies Product Datasheet

Product Name:	MYOZ2 Antibody
Catalog Number:	C31013
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human MYOZ2. AA range:21-70
Concentration:	1 mg/ml
Gene Name:	MYOZ2 C4orf5
Protein Name:	myozenin 2
Other Name:	MYOZ2 C4orf5
Human Swiss Prot No.:	Q9NPC6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	The antibody detects endogenous MYOZ2 protein
Storage and Stability:	-20°C/1 year