

Antibodies Product Datasheet

Product Name:	Plasminogen receptor Antibody
Catalog Number:	C31023
Reactivity:	Human, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human PLGRKT. AA range:1-50
Concentration:	1 mg/ml
Gene Name:	PLGRKT
Protein Name:	Plasminogen receptor
Other Name:	Plasminogen receptor (KT) (Plg-R(KT))
Human Swiss Prot No.:	Q9HBL7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Plasminogen receptor Polyclonal Antibody detects endogenous levels of Plasminogen receptor
Storage and Stability:	-20°C/1 year