

Antibodies Product Datasheet

Product Name:	NLRX1 Antibody
Catalog Number:	C31027
Reactivity:	Human, Mouse, Rat
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	WB 1:500-2000, IHC-p 1:50-300, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human NLRX1. AA range:581-630
Concentration:	1 mg/ml
Gene Name:	NLRX1
Protein Name:	NLRX1
Other Name:	NLR family member X1 (Caterpillar protein 11.3) (CLR11.3) (Nucleotide-binding oligomerization domain protein 26) (Nucleotide-binding oligomerization domain protein 5) (Nucleotide-binding oligomerization domain protein 9)
Human Swiss Prot No.:	Q86UT6
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	NLRX1 Polyclonal Antibody detects endogenous levels of NLRX1
Storage and Stability:	-20°C/1 year