

Antibodies Product Datasheet

Product Name:	CREB3 Antibody
Catalog Number:	C31028
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human CREB3. AA range:151-200
Concentration:	1 mg/ml
Gene Name:	CREB3
Protein Name:	CREB3
Other Name:	Cyclic AMP-responsive element-binding protein 3 (CREB-3) (cAMP-responsive element-binding protein 3) (Leucin zipper proitein) (Luman) (Transcription factor LZIP-alpha) [Cleaved into: Processed cyclic AMP-responsive element-binding protein 3 (N-terminal Lu
Human Swiss Prot No.:	O43889
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CREB3 Polyclonal Antibody detects endogenous levels of CREB3
Storage and Stability:	-20°C/1 year