

Antibodies Product Datasheet

Product Name:	ETBR Antibody
Catalog Number:	C31029
Reactivity:	Human, Rat, Mouse
Applications:	WB, IHC-p, IF(paraffin section), ELISA
Dilutions:	IHC-p: 100-300.WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from ETBR at AA range: 31-80
Concentration:	1 mg/ml
Gene Name:	EDNRB
Protein Name:	ETBR
Other Name:	Endothelin B receptor (ET-B) (ET-BR) (Endothelin receptor non-selective type)
Human Swiss Prot No.:	P24530
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	ETBR Polyclonal Antibody detects endogenous levels of ETBR
Storage and Stability:	-20°C/1 year