

Antibodies Product Datasheet

Product Name:	Nanos Homologue 1 (NANOS1) Antibody
Catalog Number:	C31035
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from Nanos Homologue 1 (NANOS1) at AA range: 151-200
Concentration:	1 mg/ml
Gene Name:	NANOS1
Protein Name:	Nanos Homologue 1 (NANOS1)
Other Name:	Nanos homolog 1 (NOS-1) (EC_Rep1a)
Human Swiss Prot No.:	Q8WY41
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Nanos Homologue 1 (NANOS1) Polyclonal Antibody detects endogenous levels of Nanos Homologue 1 (NANOS1)
Storage and Stability:	-20°C/1 year