

Antibodies Product Datasheet

Product Name:	Factor IX Antibody
Catalog Number:	C31036
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from Factor IX at AA range: 412-461
Concentration:	1 mg/ml
Gene Name:	F9
Protein Name:	Factor IX
Other Name:	Coagulation factor IX (EC 3.4.21.22) (Christmas factor) (Plasma thromboplastin component) (PTC) [Cleaved into: Coagulation factor IXa light chain; Coagulation factor IXa heavy chain]
Human Swiss Prot No.:	P00740
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	Factor IX Polyclonal Antibody detects endogenous levels of Factor IX
Storage and Stability:	-20°C/1 year