
Antibodies Product Datasheet

Product Name:	CLECSF6 Antibody
Catalog Number:	C31042
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the C-terminal region of human CLEC4A. AA range:161-210
Concentration:	1 mg/ml
Gene Name:	CLEC4A
Protein Name:	CLECSF6
Other Name:	C-type lectin domain family 4 member A (C-type lectin DDB27) (C-type lectin superfamily member 6) (Dendritic cell immunoreceptor) (Lectin-like immunoreceptor)
Human Swiss Prot No.:	Q9UMR7
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	CLECSF6 Polyclonal Antibody detects endogenous levels of CLECSF6
Storage and Stability:	-20°C/1 year