

Antibodies Product Datasheet

Product Name:	COL8A2 Antibody
Catalog Number:	C31055
Reactivity:	Human, Mouse, Rat
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human COL8A2. at AA range: 611-660
Concentration:	1 mg/ml
Gene Name:	COL8A2
Protein Name:	COL8A2
Other Name:	Collagen alpha-2(VIII) chain (Endothelial collagen)
Human Swiss Prot No.:	P25067
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	COL8A2 Polyclonal Antibody detects endogenous levels of COL8A2
Storage and Stability:	-20°C/1 year