

Antibodies Product Datasheet

Product Name:	C1RL Antibody
Catalog Number:	C31068
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human C1RL. at AA range: 371-420
Concentration:	1 mg/ml
Gene Name:	C1RL
Protein Name:	C1RL
Other Name:	Complement C1r subcomponent-like protein (C1r-LP) (C1r-like protein) (EC 3.4.21.-) (C1r-like serine protease analog protein) (CLSPa)
Human Swiss Prot No.:	Q9NZP8
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	C1RL Polyclonal Antibody detects endogenous levels of C1RL
Storage and Stability:	-20°C/1 year