
Antibodies Product Datasheet

Product Name:	COL17A1 Antibody
Catalog Number:	C31075
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	The antiserum was produced against synthesized peptide derived from the Internal region of human COL17A1. AA range:481-530
Concentration:	1 mg/ml
Gene Name:	COL17A1 BP180 BPAG2
Protein Name:	collagen, type XVII, alpha 1
Other Name:	COL17A1 BP180 BPAG2
Human Swiss Prot No.:	Q9UMD9
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	The antibody detects endogenous COL17A1 protein
Storage and Stability:	-20°C/1 year