

Antibodies Product Datasheet

Product Name:	NMUR2 Antibody
Catalog Number:	C31079
Reactivity:	Human, Rat, Mouse
Applications:	WB, ELISA
Dilutions:	WB 1:500-2000, ELISA 1:10000-20000
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	Synthesized peptide derived from human NMUR2. at AA range: 1-50
Concentration:	1 mg/ml
Gene Name:	NMUR2
Protein Name:	NMUR2
Other Name:	Neuromedin-U receptor 2 (NMU-R2) (G-protein coupled receptor FM-4) (G-protein coupled receptor TGR-1)
Human Swiss Prot No.:	Q9GZQ4
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.02% sodium azide.
Specificity:	NMUR2 Polyclonal Antibody detects endogenous levels of NMUR2
Storage and Stability:	-20°C/1 year