

Antibodies Product Datasheet

Product Name:	Myc-Tag Antibody, HRP Conjugated
Catalog Number:	C40004
Reactivity:	Species independent
Applications:	WB
Dilutions:	Optimal working dilutions should be determined experimentally by the investigator. Suggested starting dilutions are as follows: WB 1:5000.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit IgG
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using specific immunogen.
Immunogen:	
Concentration:	1mg/ml
Gene Name:	Myc Tag
Protein Name:	Myc-Tag
Other Name:	Myc
Human Swiss Prot No.:	
Formulation:	Liquid in PBS, pH 7.4, containing 0.02% sodium azide as preservative and 50% Glycerol.
Specificity:	Myc-Tag Polyclonal Antibody HRP Conjugated, specially designed for your Western blot analysis.
Storage and Stability:	Stable for one year at -20°C from date of shipment. For maximum recovery of product, centrifuge the original vial after thawing and prior to removing the cap. Aliquot to avoid repeated freezing and thawing.