

Antibodies Product Datasheet

Product Name:	FLAG-tag Antibody
Catalog Number:	C40009
Reactivity:	Species independent
Applications:	WB, ELISA
Dilutions:	Western Blot: 1/1000 - 1/3000. ELISA: 1/20000. Not yet tested in other applications.
Modification:	Unmodified/Total
Source:	Polyclonal Rabbit
Purification:	The antibody was affinity-purified from rabbit antiserum by affinity-chromatography using epitope-specific immunogen.
Immunogen:	DDDDK synthetic peptide conjugated to KLH.
Concentration:	1 mg/ml
Gene Name:	Flag Tag
Protein Name:	Flag-Tag DDDDK TAG, DDDDK-TAG
Other Name:	DDDDK epitope tag; DDDDK epitope tag; DYKDDDDK epitope tag
Human Swiss Prot No.:	
Formulation:	Liquid in PBS containing 50% glycerol, 0.5% BSA and 0.05% sodium azide.
Specificity:	FLAG-tag Polyclonal Antibody detects FLAG-tagged recombinant proteins or FLAG-tagged proteins overexpressed in cells.
Storage and Stability:	-20°C/1 year